


フォーカスマウント




ボール取付式フォーカスマウント

御使用の光学系や用途に合わせて様々なタイプを選択いただけます。


MS-601-AR20




MS-601-B




MS-602-AR20




MS-602-B




MS-625-295




MS-625-A




MS-625-B



MS-635-A



MS-635-B



Coarse focus mount Specification

Model	MS-601-AR20	MS-601-B	MS-602-AR20	MS-602-B		
Focusing type	Rack and Pinion type					
Slide	Cross roller system					
Stroke	±25mm					
XY stage	-	-	XY travel 13mm (Cross roller system)			
Pole diameter	Φ25mm					
Stand	MS-701,MS-703,MS-701-500					
Lens diameter	Φ29,5mm	Φ45mm	Φ29,5mm	Φ45mm		
Lens model	MS-501series,etc	MS-506,etc	MS-501series,etc	MS-506,etc		
Weight	940g	870g	1,190g	1,190g		

Coarse focus mount Specification

Model	MS-625-295	MS-625-A	MS-625-B	MS-635-A	MS-635-B
Focusing type	Rack and Pinion type				
Slide	Cross roller system				
Stroke	±15mm				
Fine adjustment	3μm/1 scale				
XY stage	Non				
Pole diameter	Φ25mm			Φ35mm	
Stand	MS-701, MS-703, MS-701-500			MS-702, 704	
Lens diameter	Φ29.5mm	-	-	-	-
Lens model	MS-501 series .etc	MS-200-TWL MS-200-TCD18 MS-200-TCD26.5 .etc	MS-200-TCD33 MS-200-TCD33-TR .etc	MS-200-TWL MS-200-TCD18 MS-200-TCD26.5 .etc	MS-200-TCD33 MS-200-TCD33-TR .etc
Weight	2,200g	2,520g	2,520g	2,700g	2,700g